1.

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

AW-CP020286 Rev E 0423

English

INDICATIONS: For sprains, strains and stable fractures of the lower leg and foot.

FITTING INSTRUCTIONS:

1. Loosen all of the straps. Detach the plastic anterior plate.

2. Open the foam liner and insert the foot and calf inside.

3. Starting with the foot section of the liner, fit the boot by overlapping the toe flap and then the cuffs. Note the toe flap is removable. One side has a hook closure that will attach over the opposite side.

4. Repeat the process for the calf section of the foam liner.

5. Re-attach the plastic anterior shell.

6. Secure the strap on the top of the foot section (1) first. Insert the strap end through the black D-ring and pull to the desired tension on the strap.

7. Next, secure the bottom calf strap (2). Insert the strap end through the black D-ring and pull to the desired tension on the strap. Repeat the process for the top calf strap (3).

7. (B) On the four strap configuration, secure bottom calf strap (2). Repeat the process for the middle (3) and top calf straps (4).

8. Push the air pump to inflate the air bladder to the desired patient comfort. The fitting process is complete.

To deflate the air bladders, press firmly on the release valve until air has been released.

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© 2023 Breg, Inc. All rights reserved.