L'ESPAGNOL

1 Desabroche la hebilla

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This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient. This device is not intended to prevent injuries or reduce or eliminate the risk of re-injury.

CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION: This device is not intended to prevent injury, or

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reduce or eliminate the risk of re-injury.

CAUTION: This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150F (65 C.) for any extended period of time. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

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








ENGLISH

1 Unfasten Buckle

- Unfasten all of the buckle ends by pushing the gray button on each buckle.
- If this is the first time the brace is applied, the brace should be at its lowest setting. If not, turn the force knob to its lowest setting.

2 Apply the Brace

Apply the brace with the patient in a seated position and the leg bent at 80°-90°.

Medial Osteoarthritis (MA): Place the hinge on the inside of the leg slightly higher and behind the joint axis (the adductor tubercle, or large bony knot on the inside of the leg).

Lateral Osteoarthritis (LA): Place the hinge on the outside of the leg slightly higher and behind the joint axis.

3 Fasten Upper Calf Strap

- Apply the upper calf strap by wrapping the strap up and over the calf muscle and around to the front of the shell to engage the buckle.
- Adjust the length of the strap using the hook to pile fastener at the buckle end of the strap.

4 Fasten Lower Calf Strap

Apply the lower calf strap in the same manner.

5 Fasten Lower Thigh Strap

- Wrap the lower thigh strap firmly around the back of the leg and engage buckle on the opposite side of the shell.
- Adjust the length of the strap using the fastener at the buckle end of the strap.

6 Fasten Upper Thigh Strap

- Wrap the upper thigh strap around and behind the leg to engage the buckle.
- Adjust the strap using the fastener at the buckle end of the strap.
- If necessary, stand and tighten straps.

7 Hinge Adjustment

Stand and walk in the brace to see how it feels on the leg. Some patients may experience pain relief prior to adjusting the hinge, due to varying leg angulations. The brace only applies force as it approaches full extension. It can only be adjusted while seated with the leg bent at 90°. If more force is required, sit down and adjust the dial to an increased force setting. Walk again to test for pain relief.

NOTE: Only adjust the hinge enough for pain relief. Walk after each adjustment. Too much force may cause discomfort and excessive muscle soreness. Once you reach pain relief, don't adjust the hinge any further. It takes at least a week for the leg to get used to having pressure applied to it. If there is any pinching from the straps, bend the leg to 90° and readjust the straps for comfort. It is important that the straps closest to the knee are snug, or you will lose brace effectiveness.

8 Locking the Hinge (Optional)

To guard against accidental movement of the force dial, use the hex wrench (provided) to lightly tighten the screw.

INDICATIONS: Thruster RLF MA is indicated for Medial Unicompartamental Osteoarthritis. Thruster RLF LA is indicated for Lateral Unicompartamental Osteoarthritis. Each brace may help to relieve symptoms of pain and inflammation. These braces are designed to aid in delaying osteotomy or total joint procedures, to increase daily walking time while increasing activities, and to aid in slowing the progression of the condition. These braces are also indicated for post-operative or post-injury application following collateral ligament injury or reconstruction (when worn on the opposite side from the injury) to relieve stress on these healing structures and post-meniscal repairs.

CONTRAINDICATIONS: Thruster RLF MA and Thruster RLF LA are contraindicated for tri- or bi-compartmental osteoarthritis, or for knees with injury or chronic stretch of the collateral ligament on the same side of the leg as the brace placement.

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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